

# Technical Data: Protecta X



Dimensions	85cm (W) x 198cm (L) x 14cm (D)
Minimum Patient Weight	No minimum patient weight
Maximum Patient Weight	250kg (39st)
Mattress Weight	14kg
Safety Standards	EN ISO 845: 2009, EN ISO 2439: 2008, EN ISO 1798: 2008, EN ISO 1856: 2018, ISO 13362:2000, EN ISO 3385:2014 & BS 3379:2005 +A1:2011
Crib Rating	The Protecta X mattress meets the requirements of BS7177:2008 for medium hazard (Crib 5)
Construction	Two foam layers. Firm base and walls with side to side profiled head and sacral zones. Geltex heel zone. Profiled base hinges. Printed X-Treme high frequency welded, moisture vapour permeable, multi-stretch cover.
Foam Grade (Base and Walls)	SB40170
Foam Grade (Head and Sacral)	SB54120
Foam Grade (Heel)	G52080
Warranty	4 years