

Technical Data:

Eleganza 5

medstrom⁺

Improved Patient Outcomes



Bed Weight	248kg
Safe Working Load	250kg (39st)
Max Patient Weight	185kg (29st)
Minimum Height	43.5cm
Maximum Height	81.5cm
Overall Bed Dimensions	100cm (W) x 219cm (L)
Sleep Deck Dimensions	90cm (W) x 208cm (L)
Bed Length Extension	22cm
Back Rest Angle	65°
Knee Rest Angle	30°
Calf Section Angle	30°
Trend/Reverse Trend Angle	+/- 14°
Lateral Tilt Angle	+/- 15°
Weigh Scales	Class III weigh scales
Castors	Four twin-wheel, 15cm multi-directional castors, with central brake and steer
CPR	Electric and bi-lateral manual
Side Rails	Split design with integrated bed controls
Warranty	24 months from purchased date

SM1148, Rev1, Feb 2023