

# Technical Data:

## Versatech 1100

**medstrom<sup>+</sup>**

Improved Patient Outcomes



Bed Weight	240kg
Safe Working Load	500kg
Minimum Height	21.6cm
Maximum Height	81.3cm
Minimum Overall Bed Dimensions	102cm (W) x 237cm (L)
Maximum Overall Bed Dimensions	132cm (W) x 257cm (L)
Sleep Deck Widths	Adjustable from 89cm/97cm/104cm/119cm
Sleep Deck Lengths	Adjustable from 203cm/213cm/223cm
Back Rest Angle	65°
Thigh Section Angle	34°
Trend/Reverse Trend Angle	+/- 12°
Under Bed Clearance	14.9cm
Castors	125mm diameter anti-static castors with steer, neutral and brake system
CPR	Electronic and manual CPR controls
Side Rails	Split side rail design

SM1178, Rev2, Jul 2023