

Technical Data:

Versatech 1100

medstrom⁺

Improved Patient Outcomes



Bed Weight (Model A)	240kg
Bed Weight (Model B)	263kg
Safe Working Load	500kg
Minimum Height	21.6cm
Maximum Height	81.3cm
Minimum Overall Bed Dimensions	102cm (W) x 237cm (L)
Maximum Overall Bed Dimensions	132cm (W) x 257cm (L)
Sleep Deck Widths	Adjustable from 89cm/97cm/104cm/119cm
Sleep Deck Lengths	Adjustable from 203cm/213cm/223cm
Back Rest Angle	70°
Thigh Section Angle	32°
Trend/Reverse Trend Angle	+/- 12°
Under Bed Clearance	14.9cm
Castors	125mm diameter anti-static castors with steer, neutral and brake system
CPR	Electronic and manual CPR controls
Side Rails	Split side rail design

SM1178, Rev3, Feb25