

# Technical Data:

## Medstrom Solo Luxe

**medstrom<sup>+</sup>**

Improved Patient Outcomes



Height Range	19cm - 81cm
Safe Working Load (SWL)	258kg (41st)
Max Patient Weight (MPW)	193kg (30st)
Back Rest Angle	70°
Knee Angle	30°
Calf Angle	20°
Trend / Reverse Trend Angle	+ 14° / - 14°
Total Length & Width	227cm x 103cm
Sleep Deck Dimensions	206cm x 91cm
Total Bed Weight	149kg (23st)

PRO00054  
CRN-0000158\_Rev1