

Technical Data:

Protecta 7



Superior profiling mattress with profiled foam core for pressure redistribution. Crib 7 fire resistant with a fully sealed and welded cover for use in mental health or very high risk environments.

Mattress dimensions	90cm (W) x 200cm (L) x 14cm (D)
Minimum patient weight	No minimum patient weight
Maximum patient weight	250kg (39st)
Mattress weight	13kg
Crib rating	BS7177:2008 for high hazard (Crib 7).
Foam specification	SB33/120. Profiled CHMR foam surface, profiled both sides.
Mattress cover	PU coated, multi-stretch, vapour permeable, waterproof with high frequency welded seams. Wipe clean cover, anti-microbial agents and resistant to bacterial growth.