

Technical Data:

AER Magna Bariatric Mattress

medstrom⁺

Improved Patient Outcomes



Mattress Dimensions (inflated)	90/105/120cm (W) x 200cm (L) x 25cm (D)
Control Unit Dimensions	34cm (W) x 24.4cm(H) x 12.5cm (D)
Maximum Patient Weight	455kg (71st)
Minimum Patient Weight	40kg
Mattress Weight	12kg
Control Unit Weight	3.5kg
Modes	Continuous low pressure, Alternating pressure, Max inflate
Comfort Setting Levels	Ten
Transport Mode	Transport cap provided, remains inflated for up to 12 hours
Alarms	Low pressure, power failure and alternate failure (visual and audible alarm)
Cover	Multi-stretch PU cover, welded seams, vapour permeable and water resistant. Concealed zipper.
Cell Material	Nylon/TPU
Warranty	2 years for unit and mattress

SM1394 Rev1.0